




## MACHEREY-NAGEL RNA Promotion

Bioanalysis



### High quality RNA purification

- Innovative products for small sample inputs
- Kits with gDNA removal columns available
- Efficient miRNA purification options
- Stabilize/preserve samples for RNA extraction

20% OFF  
plus MN hand  
sanitizer with  
your order!



**MACHEREY-NAGEL**

[www.dmarkbio.com](http://www.dmarkbio.com)

# RNA Promotion

## Featured products

### NucleoSpin® RNA Plus

High-quality RNA from cells and tissue. Efficient lysis without beta-mercaptoethanol. Filtration of lysate and gDNA removal with one column in a single step.

### NucleoZOL

Easy one-phase isolation of RNA. Separation of small (< 200 nt) and large RNA (> 200 nt) possible. Scalable protocol.

### NucleoProtect® RNA

RNA stabilization reagent for cells and tissues. Protect your samples from RNA degradation. Combine with the RNA isolation method of your choice.

### NucleoSpin® RNA XS

For use with small sample inputs. Elution volume ≥ 5 µL. rDNase included.

REF	Kit	Preps	List price (\$)	Promo Price (\$)
<b>RNA</b>				
740990.50	NucleoSpin® RNA Plus XS	50	<del>616.00</del>	492.80
740990.250	NucleoSpin® RNA Plus XS	250	<del>2,722.00</del>	2,177.60
740984.50	NucleoSpin® RNA Plus	50	<del>456.00</del>	364.80
740984.250	NucleoSpin® RNA Plus	250	<del>2,015.00</del>	1,612.00
740902.50	NucleoSpin® RNA XS	50	<del>599.00</del>	479.20
740902.250	NucleoSpin® RNA XS	250	<del>2,650.00</del>	2,120.00
740955.50	NucleoSpin® RNA	50	<del>412.00</del>	329.60
740955.250	NucleoSpin® RNA	250	<del>1,886.00</del>	1,508.80
740404.200	NucleoZOL	200 mL	<del>459.00</del>	367.20
740406.50	NucleoSpin® RNA Set for NucleoZOL	50	<del>320.00</del>	256.00
<b>miRNA</b>				
740971.50	NucleoSpin® miRNA	50	<del>546.00</del>	436.80
740971.250	NucleoSpin® miRNA	250	<del>2,386.00</del>	1,908.80
<b>RNA Stabilization</b>				
740400.50	NucleoProtect® RNA	50 mL	<del>102.00</del>	81.60
740400.250	NucleoProtect® RNA	250 mL	<del>342.00</del>	273.60
740400.500	NucleoProtect® RNA	500 mL	<del>480.00</del>	384.00

Valid until April 30, 2023. Quote RNA2023 when ordering.

20% OFF  
plus MN hand  
sanitizer with  
your order!

Promo code: RNA2023

\* Valid 2/1/2023–4/30/2023

Quote RNA2023 when you order via orders@dmartbio.com, phone, or through our Webshop.

[www.dmartbio.com](http://www.dmartbio.com)

**MACHERY-NAGEL**

DMark  
35 Grand Marshall Drive, Toronto, Ontario  
Canada M1B 5W9  
416.297.8220 | TF 1.800.784.7485



**D·Mark**

